

Analysis of Lightning and BioClean tests for assessment of sanitation levels in pork production facilities

Jason A. Kelly; Sandra F. Amass, DVM, PhD, Dipl ABVP; Darryl Ragland, DVM, PhD; Pat M. Spicer PhD; Roberta M. Alvarez, DVM

Summary

Objectives: To determine the test sensitivity and specificity of Lightning (BioControl Systems, Inc, Bellevue, Washington) and BioClean (BioVet, St Anthony, Minnesota) testing systems in determining whether various surfaces in pork production facilities have met the standards for disinfection; and to identify factors in pork production units which may interfere with the test sensitivity and specificity of Lightning or BioClean test results.

Methods: Swab samples were collected from feeders, flooring, and walls in a wean-finish room and a nursery room of

a 1600-sow, farrow-to-finish commercial operation. In both facilities, three adjacent swab samples of a 6.16 cm² (0.955 in²) area were collected at each sampling site and analyzed using Lightning, BioClean, and cultural examination for bacteria. The test sensitivity and specificity of Lightning and BioClean tests were calculated for each surface using cultural examination as the “gold standard” for classifying a sample as clean or contaminated. Factors such as feed, manure, and disinfectant residues were tested to determine if they interfered with the test sensitivity and specificity of Lightning or BioClean tests.

Results: Lightning tests were generally

highly sensitive but had low specificity. BioClean tests were generally highly specific but had low sensitivity.

Implications: Lightning and BioClean testing systems are not recommended for use in evaluating sanitation levels on swine farms without prior, independent, on-farm validation. Caution should be exercised when transferring technologies from other industries to pork production.

Keywords: swine, Lightning, BioClean, sanitation, disinfectants

Received: October 31, 2000

Accepted: May 9, 2001

All in-all out management of swine confinement units is strongly recommended to maintain herd health and increase profitability. An essential step for implementation of all in-all out pig flow is the cleaning and disinfection of rooms between pig groups to prevent exposure to pathogens shed by the previous groups.^{1,2,3,4} Recommendations include pressure washing the room and its equipment with hot water and a detergent to remove organic matter, and then disinfecting with a suitable product.^{5,6} Reports have shown that label claims do not always reflect disinfectant efficacy.⁷ Consequently, in most cases, an outbreak of neonatal scours or other costly disease must occur before producers and veterinarians reevaluate cleaning methods and disinfectants used. Currently, producers do not have an objective, rapid, proactive method to evaluate the final cleanliness of the room or the disinfectant's efficacy in inactivating pathogens enzootic to the farm.

Rapid tests to evaluate sanitation levels have not been developed specifically for use on farms. However, two testing systems, Lightning (BioControl Systems, Inc, Bellevue, Washington) and BioClean (BioVet, St Anthony, Minnesota) were developed for monitoring food safety in processing plants.^{8,9} The Lightning system has recently been used without validation for assessing cleanliness of swine transport vehicles.^{10,11}

The Lightning test method quantifies residual adenosine triphosphate (ATP) on surfaces using a luminometer and a swab sample.¹² Adenosine triphosphate is found in most food residues and bacteria. Lightning theoretically correlates increased levels of ATP with increased contamination of surfaces. The Lightning test displays a score 11 seconds after the surface sample is obtained. Each test costs approximately US \$2.78 after the initial capital outlay for the luminometer.

BioClean is a swab test that uses a color-changing reaction to quantitate protein levels on surfaces as an indicator of food residue. Increasing amounts of protein yield a more intense color change.⁸ Theoretically, surfaces with a high level of bacterial contamination have a greater protein content. The BioClean test can be prepared in approximately 1 minute and results are obtained within 20 minutes after sampling. Each BioClean test costs approximately US \$2.16.

The overall objective of the study was to determine if Lightning or BioClean tests were valid, rapid test methods to assess the efficacy of disinfection in pork production facilities. The specific objectives were to evaluate the test sensitivity and specificity of Lightning and BioClean testing systems in determining whether various surfaces in pork production facilities have met the standards for disinfection,¹ and to identify factors in pork production units which may interfere with the test sensitivity and specificity of Lightning or BioClean test results.

Our working hypothesis was that the test sensitivity and specificity of these tests would vary according to the sampled surface.

JAK, SFA, DR, RMA: Department of Veterinary Clinical Sciences, Purdue University School of Veterinary Medicine, 1248 Lynn Hall, West Lafayette, Indiana 47907-1248; PMS: Department of Natural Resources and Environment, Box 3100 Bendigo Delivery Centre, Bendigo, Victoria, Australia, 3554

This article is available online at <http://www.aasv.org/shap.html>.

Kelly JA, Amass SF, Ragland D, et al. Analysis of Lightning and BioClean tests for assessment of sanitation levels in pork production facilities. *J Swine Health Prod.* 2001;9(5):207-213.

Materials and methods

Experiment One: Determination of sensitivity and specificity of Lightning and BioClean tests on surfaces in swine facilities

Study Design

Swab samples were collected from feeders, flooring, and walls in a wean-finish room and a nursery room of a 1600-sow, farrow-to-finish commercial operation. A 26-pen wean-finish room was emptied of pigs, cleaned, disinfected with a phenolic disinfectant (Triphenol-R/256; ID Russell Company Laboratories, Longmont, Colorado), and allowed to dry for 1 day. Each pen in the facility was sampled on three surfaces: stainless steel feeder trough (Farmweld, Teutopolis, Illinois), concrete slat, and concrete wall approximately 30 cm up from the floor. Pigs were also emptied from a 16-pen nursery room. The room was cleaned and disinfected with Triphenol-R/256 and allowed to dry for 2 days. Samples were taken in each pen from the stainless steel feeder trough (Farmweld, Teutopolis, Illinois), plastic mesh flooring (MIK flooring, MIK Heinrich Michel, Marienhausen, Germany), and PVC plastic-covered wall (AP Livestock, Assumption, Illinois) approximately 30 cm up from the floor. A sterile metal washer with an internal surface area of 6.16-cm² was used to standardize the sampled surface area. In both facilities, three adjacent swab samples of a 6.16-cm² area were collected at each sampling site and analyzed using Lightning, BioClean, and cultural examination for aerobic bacteria.

Background contamination control samples

Sterile washers were exposed to facility airspace for approximately 5 seconds at five different locations in each room to quantify background aerosol contamination of washers that might occur during sampling. The inner surface of each washer was swabbed and analyzed using Lightning, BioClean, and cultural examination for bacteria as described below. Additionally, metal washers autoclaved in the same batch as those used for on-farm sampling were culturally examined under aseptic conditions to confirm sterility.

Bacterial plating analysis

Each surface was sampled with sterile, cotton-tipped swabs which were placed into 1-mL aliquots of 0.9% saline stored on ice

Table 1: Calculation of test parameters for two testing systems, Lightning¹ and BioClean², used to assess effectiveness of cleaning and disinfecting surfaces in swine facilities. Cultural examination was performed concurrently on the same surfaces for each test.

Test parameter	Calculation
Test sensitivity	$\frac{\text{No. true positive test results determined by the testing system}^3}{\text{No. positives determined by culture}}$
Test specificity	$\frac{\text{No. true negative test results determined by the testing system}^3}{\text{No. negatives determined by culture}}$
Test positive predictive value	$\frac{\text{No. true positive test results determined by the testing system}}{\text{Total no. positive results determined by the testing system}^4}$
Test negative predictive value	$\frac{\text{No. true negative results determined by the testing system}}{\text{Total no. negative results determined by the testing system}^4}$

¹ BioControl Systems, Inc, Bellevue, Washington

² BioVet, St Anthony, Minnesota

³ True positive test results are positive (contaminated) by culture; true negative test results are negative (clean) by culture.

⁴ Total no. of positive test results includes results that are positive (contaminated) by the testing system and either positive or negative (clean) by culture; total no. of negative test results includes results that are negative by the testing system and either positive or negative by culture.

packs. Original samples and serial dilutions were plated onto trypticase soy agar with 5% sheep blood (BBL Stacker Plate, Beckton Dickenson Microbiology Systems, Cockeysville, Maryland) within 5 hours of collection and incubated aerobically at 37°C for 18 hours. Colonies were counted on each plate and converted to colony forming units (CFUs) per cm² of surface sampled. The average bacterial count for the five background contamination control samples for the room was subtracted from each sample swab count. A surface was classified as either clean (≤ 1 CFU per cm²) or contaminated (> 1 CFU per cm²).¹ Plates indicating more than 487 CFUs per cm² surface area were considered “too numerous to count” (TNTC) and classified as contaminated. Samples that detected less than the mean background contamination level were classified as sterile (0 CFUs per cm²).

Lightning analysis

The Lightning system was used according to test instructions. Briefly, each surface was sampled with a Lightning swab that consists of three basic parts. A swab, attached to an ampule containing buffer solution at the top, is enclosed in a sheath containing a luciferin-luciferase pellet. Buffer contained in the ampule is used to

rinse the sample from the swab and mix with the luciferin-luciferase pellet to activate the reaction. Residual ATP present on surfaces reacts with the luciferin-luciferase pellet to yield light. The Lightning luminometer measures the light output and calculates a score. Each sample was classified by the standard Lightning cut-off as either clean (score ≤ 2.5) or contaminated (score > 2.5). Alternative cut-off values (3.0, 3.5, 4.0) were applied to determine their effects on test sensitivity and specificity.

BioClean analysis

Sample tubes were prepared according to test instructions by adding 1 drop of reagent B to the tube containing reagent A. Each surface was sampled with a sterile cotton swab provided with the BioClean kit. Protein contaminants react with copper ions to form a complex with the biuret reagent that causes a color-changing reaction. The sample swab was allowed to react for 20 minutes in the tube containing the reagent mixture. The surface was classified as clean (no color change after 20 minutes) or contaminated (any color change after 20 minutes).

Data analysis

The sensitivity, specificity, positive predic-

tive value (PPV), and negative predictive value (NPV) of Lightning and BioClean tests were calculated (Table 1) for each surface, using cultural examination as the “gold standard” for classifying a sample as truly clean or contaminated. Sensitivity was defined as the ability of the test to detect a contaminated surface when the surface was truly contaminated. Specificity was defined as the ability of the test to detect a clean surface when the surface was truly clean. The PPV was defined as the probability of a contaminated test result corresponding to a truly contaminated surface. The NPV was defined as the probability of a clean test result corresponding to a truly clean surface.

Experiment Two: Factors in pork production units that may interfere with sensitivity and specificity of Lightning or BioClean tests

Effect of manure residue on test results

Five manure samples were collected from a growing-finishing facility housing pigs weighing approximately 45 kg. Manure samples were diluted with sterile 0.22- μ m-filtered water to a final dilution of 0.001 g manure per mL water. Each sample was

divided into two equal aliquots. One aliquot was autoclaved for 20 minutes at 121°C and the other remained nonautoclaved. A negative control consisted of a sterile, 0.22- μ m-filtered water sample. The nonautoclaved manure samples, the autoclaved manure samples, and a negative control were analyzed using cultural examination for bacteria, Lightning, and BioClean methods.

Effect of feed residue on test results

Five grow-finish feed samples were collected. Feed samples were diluted with sterile 0.22- μ m-filtered water to a final dilution of 0.083 g feed per mL water. Each sample was divided into two equal aliquots. One aliquot was autoclaved for 20 minutes at 121°C and the other remained nonautoclaved. A negative control consisted of a sterile, 0.22- μ m-filtered water sample. The nonautoclaved feed samples, the autoclaved feed samples, and a negative control were analyzed using cultural examination for bacteria, Lightning, and BioClean tests.

Effect of disinfectant residue on test results

Disinfectant sampling was performed

under controlled laboratory conditions to determine if disinfectant residue affected Lightning or BioClean test results. Representatives of six classes of disinfectants were used (Table 2). Clean glass slides (Esco 3”x 1” microscope slides, Erie Scientific Company, Portsmouth, New Hampshire) were submerged in disinfectant solution prepared according to label instructions in sterile, 0.22- μ m-filtered water. Slides were dried in a sterile, HEPA-filtered biosafety cabinet. Dry slides were sampled using each of the three testing methods. The top one third of the slide was culturally examined for bacteria, the middle one third was sampled for the Lightning test, and the bottom one third was sampled for the BioClean test. Five replicates were performed for each disinfectant. Five control slides were submerged in sterile, 0.22- μ m-filtered water and tested as above.

Effects of disinfected manure and feed on test results

Five manure samples were collected from a continuous-flow finishing facility housing pigs weighing approximately 68 kg. Manure samples were diluted in sterile, 0.22- μ m-filtered water to a concentration of 0.01 g manure per mL water. Five samples

Table 2: Active ingredients and dilution rates of disinfectants applied to clean surfaces in swine facilities to determine the test sensitivities and specificities of two testing systems, Lightning¹ and BioClean², for assessing effectiveness of disinfecting procedures.

Disinfectant name and manufacturer	Category	Active ingredient(s)	Dilution
Clorox Bleach The Clorox Company Oakland, California	chlorine	sodium hypochlorite	1 cup/gal (46.9 mL/L)
1-Stroke Environ STERIS Corporation St. Louis, Missouri	phenol	o-phenylphenol o-benzyl-p-chlorophenol p-tertiary-amyphenol	0.5oz/gal (3.9 mL/L)
Nolvasan solution Fort Dodge Laboratories Fort Dodge, Iowa	chlorhexidine	1-1' -hexamethylenebis [5-(p-chlorophenyl) biguanide] diacetate	1oz/gal (7.8 mL/L)
Roccal-D Plus Pharmacia and Upjohn Kalamazoo, Michigan	quaternary ammonium	didecyl dimethyl ammonium chloride alkyl (C ₁₂ , 61%; C ₁₄ , 23%; C ₁₆ , 11%; C ₁₈ , 2.5%; C ₈ & C ₁₀ , 2.5%) dimethyl benzyl ammonium chloride alkyl (C ₁₂ , 40%; C ₁₄ , 50%; C ₁₆ , 10%; C ₁₈ , 10%) dimethyl benzyl ammonium chloride bis-n-tributyltin oxide	0.5oz/gal (3.9 mL/L)
Cidex Activated Dialdehyde solution Advanced Sterilization Products Irvine, California	aldehyde	glutaraldehyde	Undiluted
VirkonS Farnam Livestock Products Phoenix, Arizona	oxidizing agent	potassium peroxymonosulfate sodium chloride	34.1g/gal (9.0 g/L)

¹ BioControl Systems, Inc, Bellevue, Washington

² BioVet, St Anthony, Minnesota

Table 3: Classification of surfaces as clean¹ or contaminated² by each of three test methods performed in triplicate on clean, disinfected surfaces in swine facilities.

Surface	n	True value (Cultural examination for bacteria)		Lightning ³ (cut-off score = 2.5)		BioClean ⁴	
				Contaminated	Clean	Contaminated	Clean
Stainless steel wean-finish feeders	26	Contaminated	20	20	0	3	17
		Clean	6	6	0	0	6
Stainless steel nursery feeders	16	Contaminated	10	8	2	5	5
		Clean	6	3	3	1	5
Concrete slats	26	Contaminated	18	17	1	3	15
		Clean	8	8	0	1	7
Plastic flooring	16	Contaminated	11	10	1	4	7
		Clean	5	4	1	5	0
Concrete walls	26	Contaminated	18	18	0	2	16
		Clean	8	8	0	0	8
PVC plastic-covered walls	16	Contaminated	1	0	1	0	1
		Clean	15	3	12	3	12

¹ Surfaces were classified as clean (negative test) if the average count was < 1 colony forming unit (CFU) per cm²

² Surfaces were classified as contaminated (positive test) if the average bacterial count was ≥1 CFU per cm²

³ BioControl Systems, Inc, Bellevue, Washington

⁴ BioVet, St Anthony, Minnesota

Table 4: Sensitivity, specificity, positive predictive value (PPV) and negative predictive value (NPV) of Lightning¹ (cut-off score = 2.5) and BioClean² tests on various clean, disinfected³ surfaces in a swine facility.

Surface	Test	Sensitivity (%)	Specificity (%)	PPV (%)	NPV (%)
Stainless steel wean-finish feeders	Lightning	100	0	77	NA ⁴
	BioClean	15	100	100	26
Stainless steel nursery feeders	Lightning	80	50	73	60
	BioClean	50	83	83	50
Concrete slats	Lightning	94	0	68	0
	BioClean	17	88	75	32
Plastic flooring	Lightning	91	20	71	50
	BioClean	36	0	44	0
Concrete walls	Lightning	100	0	69	NA ⁴
	BioClean	11	100	100	33
PVC plastic-covered walls	Lightning	0	80	0	92
	BioClean	0	80	0	92

¹ BioControl Systems, Inc, Bellevue, Washington

² BioVet, St Anthony, Minnesota

³ Clean surfaces disinfected with a phenolic compound and allowed to dry for 1 to 2 days prior to testing.

⁴ Not applicable: negative predictive values cannot be calculated when a test does not provide negative results

Table 5: Sensitivity and specificity of Lightning¹ test results at varied cut-off points when tests were performed on clean, disinfected² surfaces in swine facilities.

Surface	Cut-off score = 2.5		Cut-off score = 3.0		Cut-off score = 3.5		Cut-off score = 4.0	
	Sensitivity (%)	Specificity (%)	Sensitivity (%)	Specificity (%)	Sensitivity (%)	Specificity (%)	Sensitivity (%)	Specificity (%)
Stainless steel wean-finish feeders (n=26)	100	0	80	0	40	33	25	67
Stainless steel nursery feeders (n=16)	80	50	70	50	60	50	20	50
Concrete slats (n=26)	94	0	94	13	61	50	17	100
Plastic flooring (n=16)	91	20	82	60	55	60	27	100
Concrete walls (n=26)	100	0	78	0	17	63	6	100
PVC plastic-covered walls (n=16)	0	80	0	100	0	100	0	100

¹ BioControl Systems, Inc, Bellevue, Washington

² Clean surfaces disinfected with phenolic compound and allowed to dry for 1 to 2 days prior to testing.

of an SEW diet containing carbadox were collected and diluted in sterile, 0.22- μ m-filtered water to a concentration of 0.01 g feed per mL water. To simulate contamination with organic matter and bacteria, seven slides were submerged in each sample of manure and seven slides were submerged in each sample of feed. All slides were allowed to dry in a sterile, HEPA-filtered biosafety cabinet. Then, for both the manure-coated slides and the feed-coated slides, six of the seven inoculated slides were each submerged in one of the six disinfectants (Table 2) for 5 seconds and dried a second time. The seventh slide was used as a positive control that was dried but not disinfected. Five replicates were performed for each disinfectant and positive control. A negative control slide was submerged in sterile 0.22- μ m-filtered water and allowed to dry in the biosafety cabinet to confirm sterility of the water and the biosafety cabinet. All slides were evaluated using cultural examination for bacteria, Lightning, and BioClean tests. The top one third of the slide was culturally examined for bacteria, the middle one third was sampled for the Lightning test, and the bottom one third was sampled for the BioClean test.

Results

Experiment One: Determination of sensitivity and specificity of Lightning and BioClean tests on surfaces in swine facilities

Variation existed among cultural examination for bacteria, Lightning, and BioClean tests with respect to classification of surfaces as clean or contaminated (Table 3).

Lightning was highly sensitive but had low specificity on stainless steel wean-finish feeders, stainless steel nursery feeders, concrete slats, plastic flooring, and concrete walls. However, sensitivity was low and specificity was high on PVC plastic-covered walls using the standard cut-off score of 2.5 to distinguish between clean and contaminated surfaces (Table 4). Incrementally increasing the Lightning cut-off score did not improve the accuracy of test results on surfaces (Table 5). BioClean had high specificity and low sensitivity on stainless steel wean-finish feeders, stainless steel nursery feeders, concrete slats, concrete walls, and PVC plastic-covered walls. BioClean sensitivity on plastic flooring was low but exceeded specificity (Table 4).

Experiment Two: Factors in pork production units that may interfere with sensitivity and specificity of Lightning or BioClean tests

Effect of manure residue on test results

Negative control: A sterile water negative control was classified as clean using all three tests.

Nonautoclaved manure: All nonautoclaved manure samples grew bacteria in numbers TNTC on cultural examination and were classified as contaminated. Lightning (cut-off score 2.5) correctly classified four of five samples (80%) as contaminated and incorrectly classified one of five samples (20%) as clean. BioClean incorrectly classified five of five samples (100%) as clean.

Autoclaved manure: All autoclaved manure samples were sterile on cultural examina-

tion and classified as clean. Lightning (cut-off score 2.5) correctly classified three of five samples (60%) as clean and incorrectly classified two of five samples (40%) as contaminated. BioClean correctly classified all autoclaved manure samples as clean.

Effect of feed residue on test results

Negative control: A sterile water negative control was classified as clean using all three tests.

Nonautoclaved feed: All nonautoclaved feed samples grew bacteria in numbers TNTC on cultural examination and were classified as contaminated. Both Lightning (cut-off score 2.5) and BioClean tests correctly classified all samples as contaminated.

Autoclaved feed: All autoclaved feed samples were sterile on cultural examination and classified as clean. Both Lightning (cut-off score 2.5) and BioClean tests incorrectly classified five of five autoclaved feed samples (100%) as contaminated.

Effect of disinfectant residue on test results

Thirty-five of 35 slides (100%) were sterile on cultural examination. Lightning (cut-off score 2.5) correctly classified all disinfected slides and control slides as clean. BioClean correctly classified all control slides and all slides coated with Clorox Bleach, Nolvasan solution, 1-Stroke Environ, Roccal D-Plus, and Virkon S as clean. However, BioClean incorrectly classified all slides coated with Cidex Activated Dialdehyde solution as contaminated.

Effect of disinfected manure on test results

Negative control: A negative control slide coated with sterile water was classified as clean using all three tests.

Disinfected manure: Twenty-five of 25 disinfected manure samples (100%) were sterile on cultural examination and classified as clean. Lightning (cut-off score 2.5) correctly classified as clean three of five (60%) manure-coated slides disinfected with Cidex Activated Dialdehyde solution and incorrectly classified two of five (40%) as contaminated. Lightning (cut-off score 2.5) correctly classified one of five (20%) manure-coated slides disinfected with Clorox Bleach as clean and incorrectly classified four of five (80%) as contaminated. Lightning (cut-off score 2.5) correctly classified four of five (80%) manure-coated slides disinfected with Nolvasan solution as clean and incorrectly classified one of five (20%) as contaminated. Lightning (cut-off score 2.5) correctly classified one of five (20%) manure-coated slides disinfected with 1-Stroke Environ as clean and incorrectly classified four of five (80%) as contaminated. Lightning (cut-off score 2.5) correctly classified two of five (40%) manure-coated slides disinfected with Roccal D-Plus as clean and incorrectly classified three of five (60%) as contaminated. BioClean correctly classified all disinfected manure samples as clean.

Effect of disinfected feed on test results

Negative controls: A control slide coated with sterile water was classified as clean using all three tests.

Disinfected feed: Twenty-five of 25 disinfected feed samples (100%) were sterile on cultural examination and classified as clean. Lightning (cut-off score 2.5) correctly classified all feed-coated slides disinfected with Cidex Activated Dialdehyde solution or Clorox Bleach as clean. Lightning (cut-off score 2.5) correctly classified two of five (40%) feed-coated slides disinfected with Nolvasan solution as clean and incorrectly classified three of five (60%) as contaminated. Lightning (cut-off score 2.5) correctly classified three of five (60%) of feed-coated slides disinfected with 1-Stroke Environ or Roccal D-Plus as clean and incorrectly classified two of five (40%) as contaminated. BioClean correctly classified all disinfected feed samples as clean.

Discussion

An ideal test for assessing sanitation on pork production facilities would be inexpensive, rapid, 100% sensitive, and 100% specific. Lightning and BioClean testing systems were not designed for use in pork production facilities. Lightning and BioClean technology were transferred to swine production because of the need by the industry to quickly and objectively evaluate sanitation protocols. Unfortunately, Lightning tests were recommended and implemented without prior validation for on-farm use.^{10,11}

To compensate for the greater contamination of swine facilities compared to food processing plants, some have recommended using a higher Lightning cut-off score for classifying a surface as clean. In our study, incremental increases in cut-off scores improved specificity, but decreased sensitivity of Lightning for all surfaces except PVC plastic-covered walls. Therefore, altering cut-off scores did not improve overall test accuracy.

Further evaluation of Lightning under experimental conditions demonstrated that sterile organic material such as autoclaved or disinfected feed and manure residues were sometimes classified as contaminated according to the Lightning tests. Feed is composed primarily of plant products. Living cells present in the ground seed coats may contain sufficient ATP to cause false-positive Lightning results. Sterile manure residue could cause false-positive results because of ATP in fibrous seed coats, which pass through the digestive tract of the pig with little degradation. An additional source of ATP in manure could be epithelial cells sloughed from the digestive tract. Pure disinfectant residues did not affect Lightning results.

BioClean tests were recently made available for commercial use and have not been widely implemented in the pork industry. Under the conditions of this study, BioClean tests were generally highly specific for all surfaces except plastic flooring, but sensitivity was low. Both sensitivity and specificity were low for plastic flooring. Low overall BioClean sensitivity might result from inability of reagents in the kit to adequately detect small quantities of protein. An alternative explanation is that residual protein is not a good indicator of bacterial contamination. A second disad-

vantage of the BioClean test is cross-reaction with Cidex Activated Dialdehyde solution, an aldehyde disinfectant, resulting in the classification of sterile solutions of pure disinfectant as contaminated. However, false classification did not occur when Cidex Activated Dialdehyde solution was used to disinfect feed or manure samples. One explanation is that the organic material bound to or reacted with the residue in Cidex that caused false-positive reactions. Other classes of disinfectants tested did not affect BioClean test results.

In conclusion, the low test specificity of Lightning and the low test sensitivity of BioClean tests resulted in an inability to accurately assess the efficacy of disinfection in swine production facilities under the conditions of this study. This experiment took place in two rooms of a single commercial pork production facility. The walls, equipment, and flooring of this facility may not represent those on other swine farms. The investigators encourage practitioners and producers to validate Lightning and BioClean testing systems in their own facilities before implementing biosecurity programs utilizing these testing methods. Future research should focus on validating additional rapid testing systems used by other industries to determine whether they are applicable to the pork industry. Development of rapid testing systems designed for use in the pork industry should be considered if technology transfer from other industries is not possible.

Implications

- Lightning and BioClean testing systems are not recommended for use in evaluating sanitation on swine farms unless prior, independent, on-farm validation is performed.
- Caution should be exercised when transferring technologies from other industries to pork production.
- When evaluated on three different surface types in a wean-finish building and a nursery, under the conditions of this study, Lightning tests were generally highly sensitive but specificity was low when used on stainless steel feeders, concrete slats, and plastic flooring, and Lightning tests were highly specific but had low sensitivity when used on PVC plastic-covered walls.

- When evaluated on three different surface types in a wean-finish building and a nursery, under the conditions of this study, BioClean tests were highly specific, but sensitivity was low.
- Under the conditions of this study, Lightning tests classified sterile organic material such as autoclaved or disinfected feed and manure residues as contaminated, and BioClean tests cross-reacted with an aldehyde disinfectant resulting in false classification of sterile disinfectant solutions as contaminated.

Acknowledgements

Support for this research was provided by the National Pork Producers Council.

References — refereed

1. Tamasi G. Testing disinfectants for efficacy. *Rev Sci Tech.* 1995;14:75–79.
2. Quinn PJ. Disinfection and disease prevention in veterinary medicine. In: Block, SS. *Disinfection, sterilization, and preservation.* 4th ed. Philadelphia: Lea and Febiger; 1991:846–868.
3. Fotheringham VJC. Disinfection of livestock production premises. *Rev Sci Tech.* 1995;14:191–205.
4. Owen JM. Disinfection of farrowing pens. *Rev Sci Tech.* 1995;14:381–391.
5. Bruins G, Dyer JA. Environmental considerations of disinfectants used in agriculture. *Rev Sci Tech.* 1995;14:81–94.
7. Kennedy MA, Mellon VS, Caldwell G, Potgieter LND. Virucidal efficacy of the newer quaternary ammonium compounds. *JAVMA.* 1995;31:254–257.
12. Ford SR, Chenault KD, Hall MS, Pangburn SJ, Leach FR. Effect of periodate-oxidized ATP and

other nucleotides on firefly luciferase. *Arch Biochem Phys.* 1994;314:261–267.

References — non refereed

6. Thompson, R. Transportation cleaning and disinfection. *Swine Health Fact Sheet.* Jan. 2000; 2(2).
8. BioClean qualitative control of surface cleanliness: Your surfaces are clean or Bioclean? (Biovet, St Anthony, Minnesota)
9. Lightning Luminometer Operator's Manual. Westbrook, Maine: IDEXX Laboratories, Inc; 35.
10. Faust, C. The determination of surface levels of ATP as a biosecurity measure. *Proc AASP.* Quebec City, Quebec. 1997;301–311.
11. Dee, SA. An overview of methods for measuring the impact of sanitation procedures for swine transport vehicles. *Swine Health Fact Sheet.* November 1999;1(2).

